

Anti-Nestin / NES antibody



Catalog Number: 100855

Product name

Anti-Nestin / NES antibody

Specificity

Human Nestin

Antibody description

Rabbit polyclonal to Nestin / NES

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the C-terminus of the Human Nestin, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC-P

Dilutions


IHC-P: 0.1-2 µg/mL

Validations



Nestin Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of human Nestin in human brain with rabbit polyclonal antibody (1 µg/mL, formalin-fixed paraffin embedded



sections).

Nestin Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of human Nestin in human malignant melanoma (from 4 donors) with rabbit polyclonal antibody (1 µg/mL, formalin-fixed paraffin embedded sections).