

Anti-APOOL antibody



Catalog Number: 175855

Product name

Polyclonal

Anti-APOOL antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat

Applications

Antibody description

WB

Rabbit polyclonal antibody to APOOL

Dilutions

Preparation

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human apool. aa range:219-268

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Validations




Fig1.; Western blot analysis of SH-SY5Y 293T 3T3 lysis using ApoOL antibody. Antibody was diluted at 1:1000