

AIM2 antibody

Catalog Number: 107855

Product name

AIM2 antibody

Specificity

Human; other species not tested.

Antibody description

AIM2 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, A549 cells, COLO 320 cells, HeLa cells, HepG2 cells, Jurkat cells, MCF7 cells, SH-SY5Y cells. Positive IHC detected in human testis tissue. Observed molecular weight by Western-blot: 40-45kd

Preparation

This antibody was obtained by immunization of AIM2 recombinant protein (Accession Number: BC010940). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:107855(AIM2 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human testis using Catalog No:107855(AIM2 antibody) at dilution of 1:100 (under 10x lens)