

# GRAMD1B antibody



Catalog Number: 111155

## Product name

GRAMD1B antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

GRAMD1B Rabbit Polyclonal antibody. Positive WB detected in A549 cells. Observed molecular weight by Western-blot: 85 kDa

## Preparation

This antibody was obtained by immunization of GRAMD1B recombinant protein (Accession Number: NM\_020716). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:111155(GRAMD1B Antibody) at dilution of 1:1000