

Apoe Polyclonal Antibody

Catalog Number: 163155

Product name

Apoe Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to Apoe

Preparation

Antigen: A synthetic peptide of mouse Apoe

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:200

Validations



Western blot - Apoe Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using Apoe antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunofluorescence - Apoe Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using Apoe antibody at dilution of 1:100. Blue: DAPI for nuclear staining.