

Human RAB27B/C25KG (His Tag) recombinant protein



Catalog Number: 503155

General Information

Gene Name Synonym

C25KG

Protein Construction

A DNA sequence encoding the human RAB27B (NP_004154.2) (Thr2-Lys215) was expressed with an N-terminal polyhistidine tag.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The recombinant human RAB27B comprises 234 amino acids and has a predicted molecular mass of 26.6 kDa. The apparent molecular mass of the protein is approximately 29 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human RAB27B Protein (His Tag) SDS-PAGE