

# NFKB1,p105,p50-Specific antibody



Catalog Number: 113155

---

## Product name

NFKB1,p105,p50-Specific antibody

## Specificity

Human; other species not tested.

## Antibody description

NFKB1,p105,p50-Specific Rabbit Polyclonal antibody. Positive WB detected in A431 cells, Jurkat cells, Raji cells. Observed molecular weight by Western-blot: 50kd,105kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_003998). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:113155(NFKB1,p105,p50-Specific antibody) at dilution of 1:200