

CHOP; GADD153 antibody



Catalog Number: 107155

Product name

CHOP; GADD153 antibody

Specificity

Human; other species not tested.

Antibody description

CHOP; GADD153 Mouse Monoclonal antibody. Positive WB detected in Tunicamycin treated HeLa cells. Observed molecular weight by Western-blot: 30 kDa

Preparation

This antibody was obtained by immunization of CHOP; GADD153 recombinant protein (Accession Number: NM_004083). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Tunicamycin treated HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107155(CHOP; GADD153 Antibody) at dilution of 1:1000