

Anti-Cytokeratin 18 (2D6) Mouse antibody



Catalog Number: 168155

Product name

WB

Anti-Cytokeratin 18 (2D6) Mouse antibody

Dilutions

Specificity

WB: 1/1000-3000

Human, Mouse, Rat

Antibody description

Cytokeratin 18 (2D6) Mouse monoclonal antibody

Preparation

Antigen: Recombinant Protein of Keratin type I cytoskeletal 18

Formulation

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality


Monoclonal

Ig Type


Mouse IgG1

Applications

Validations



Western blot detection of Keratin 18 in HeLa, 3T3 cell lysates using Keratin 18 (2D6) Mouse mAb (1:1000 diluted). Predicted band size: 46 kDa. Observed band size: 46 kDa.



Western blot analysis of 1) HepG2, 2) HeLa, 3) Mouse Liver tissue, 4) C2C12, 5) Rat Heart tissue, 6) Mouse Skeletal Muscle tissue, (8F2) diluted at 1:2000.