

Human 14-3-3 beta/YWHAB (GST Tag) recombinant protein



Catalog Number: 502555

General Information

Gene Name Synonym

Protein 1054; Protein kinase C inhibitor protein 1; 14-3-3 protein beta/alpha, N-terminally processed

Protein Construction

A DNA sequence encoding the human YWHAB (NP_003395.1) (Met 1-Asn 246) was fused with the GST tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

> 96 % as determined by SDS-PAGE

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant human YWHAB/GST chimera

consists of 478 amino acids and has a predicted molecular mass of 54.9 kDa. It migrates as an approximately 52 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile 20mM Tris, 150mM NaCl, 0.1mM DTT, 10% glycerol, pH 7.5

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human 14-3-3 beta / YWHAB Protein (GST Tag)
SDS-PAGE