

# SPATA7 antibody

Catalog Number: 115555

## Product name

SPATA7 antibody

IP: 1:200-1:1000

IHC: 1:20-1:200

## Specificity

Human, Mouse; other species not tested.

## Antibody description

SPATA7 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, PC-3 cells, Y79 cells. Positive IP detected in mouse testis tissue. Positive IHC detected in human testis tissue, human gliomas tissue. Observed molecular weight by Western-blot: 68-80 kDa

## Preparation

This antibody was obtained by immunization of SPATA7 recombinant protein (Accession Number: BC008656). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IP, IHC

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000


## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:115555(SPATA7 antibody) at dilution of 1:300



IP Result of anti-SPATA7 (IP:Catalog No:115555, 3ug; Detection:Catalog No:115555 1:300) with mouse testis tissue lysate 8000ug.




Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog

# SPATA7 antibody

Catalog Number: 115555

---

No:115555(SPATA7 Antibody) at dilution of 1:50  
(under 10x lens)



Immunohistochemistry of paraffin-embedded  
human testis tissue slide using Catalog  
No:115555(SPATA7 Antibody) at dilution of 1:50  
(under 40x lens)