

# Anti-Transferrin Receptor Rabbit antibody



Catalog Number: 168555

## Product name

Anti-Transferrin Receptor Rabbit antibody

## Ig Type

Rabbit IgG

## Specificity

Human

## Applications

WB, IP

## Antibody description

Transferrin Receptor Rabbit monoclonal antibody

## Dilutions

1:2000

## Preparation

1:20

Antigen: A synthetic peptide of human Transferrin Receptor

## Validations

## Formulation


Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal



Western blot detection of Transferrin Receptor in HeLa cell lysates using Transferrin Receptor Rabbit mAb(1:1000 diluted). Predicted band size:84kDa. Observed band size:90kDa.