

# HAPLN1 Polyclonal Antibody



Catalog Number: 164955

## Product name

HAPLN1 Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to HAPLN1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 165-354 of human HAPLN1 (NP\_001875.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:10 - 1:100

## Validations



## Western blot - HAPLN1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using HAPLN1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



## Immunohistochemistry - HAPLN1 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human liver cancer using HAPLN1 antibody at dilution of 1:100 (40x lens).

# HAPLN1 Polyclonal Antibody




Catalog Number: 164955

---



Immunohistochemistry - HAPLN1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using HAPLN1 antibody at



dilution of 1:100 (40x lens).

Immunofluorescence - HAPLN1 Polyclonal Antibody

Immunofluorescence analysis of A549 cells using HAPLN1 antibody . Blue: DAPI for nuclear staining.