

Anti-Annexin VI Rabbit antibody

Catalog Number: 166955

Product name

Anti-Annexin VI Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Annexin VI Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human Annexin VI

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, ICC/IF

Dilutions


WB: 1/1000

ICC/IF: 1/20

Validations



Immunocytochemistry of Annexin VI(green) in Hela cells using Annexin VI Rabbit pAb at dilution



Western blot detection of Annexin VI in K562, Rat Brain, C6, 3T3 cell lysates using Annexin VI Rabbit pAb(1:1000 diluted). Predicted band size: 76kDa. Observed band size: 70kDa.