

# NR1D1 Polyclonal Antibody



Catalog Number: 162255

## Product name

Rabbit IgG

NR1D1 Polyclonal Antibody

## Applications

## Specificity

WB, IP

Human, Mouse

## Dilutions

## Antibody description

WB 1:500 - 1:2000

IP 1:20 - 1:50

Polyclonal antibody to NR1D1

## Preparation

Antigen: Recombinant protein of human NR1D1

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage


Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

## Validations



Western blot - NR1D1 Polyclonal Antibody

Western blot analysis of extracts of A-549 cells, using NR1D1 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.