

# PSTPIP1 antibody



Catalog Number: 114255

## Product name

PSTPIP1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PSTPIP1 Rabbit Polyclonal antibody. Positive WB detected in mouse thymus tissue. Positive IHC detected in human brain tissue. Observed molecular weight by Western-blot: 45-48kd

## Preparation

This antibody was obtained by immunization of PSTPIP1 recombinant protein (Accession Number: NM\_003978). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

## Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:114255(PSTPIP1 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human brain using Catalog No:114255(PSTPIP1 antibody) at dilution of 1:50 (under 40x lens)