

ZNRF2-Specific antibody



Catalog Number: 117255

Product name

ZNRF2-Specific antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ZNRF2-Specific Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, mouse brain tissue, mouse testis tissue. Observed molecular weight by Western-blot: 30-35 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_147128). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:117255(ZNRF2-Specific Antibody) at dilution of 1:300