

IL10RA antibody



Catalog Number: 111655

Product name

IL10RA antibody

Specificity

Human; other species not tested.

Antibody description

IL10RA Rabbit Polyclonal antibody. Positive WB detected in human liver tissue, HeLa cells, human heart tissue, human placenta tissue, K-562 cells. Positive FC detected in K-562 cells. Observed molecular weight by Western-blot: 70-72kd,90-95kd

Preparation

This antibody was obtained by immunization of IL10RA recombinant protein (Accession Number: NM_001558). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, FC

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:111655(IL10RA antibody) at dilution of 1:400 Plot P03, gated on P01.R1



1X10⁶ K-562 cells were stained with 0.2ug IL10RA antibody (Catalog No:111655, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1000.