

Anti-Lrp1 antibody

Catalog Number: 176655

Product name

Anti-Lrp1 antibody

Specificity

Mouse

Antibody description

Mouse monoclonal antibody to Lrp1

Preparation

This antigen of this antibody was recombinant protein

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications

WB, IHC-P

Dilutions

WB: 1:1,000

IHC-P: 1:50

Validations





Fig1: Western blot analysis on recombinant protein using anti-LRP-1 Mouse mAb (Cat. #



176655#).

Fig2: Immunohistochemical analysis of paraffin-embedded mouse liver tissue using anti-LRP-1 Mouse mAb (Cat. # 176655#).