

Anti-HNF-4-alpha Rabbit antibody



Catalog Number: 168655

Product name

Rabbit IgG

Anti-HNF-4-alpha Rabbit antibody

Applications

Specificity

WB, IP, ChIP

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1/1000

HNF-4-alpha Rabbit monoclonal antibody

IP: 1/20

Preparation


ChIP: 1/20

Antigen: Recombinant protein of human HNF-4-alpha

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Western blot detection of HNF-4-alpha in K562,C6,HeLa cell lysates using HNF-4-alpha Rabbit mAb(1:1000 diluted).Predicted band size:53kDa.Observed band size:53kDa.

Ig Type