

# NUDT11 antibody



Catalog Number: 113355

## Product name

NUDT11 antibody

## Specificity

human, mouse; other species not tested.

## Antibody description

NUDT11 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, human brain tissue, mouse testis tissue. Positive IHC detected in human ovary tumor. Observed molecular weight by Western-blot: 19 kDa

## Preparation

This antibody was obtained by immunization of NUDT11 recombinant protein (Accession Number: NM\_018159). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113355(NUDT11 antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:113355(NUDT11 antibody) at dilution of 1:100 (under 10x lens)