

# Influenza B HA (His Tag) recombinant protein



Catalog Number: 504355

## General Information

### Protein Construction

A DNA sequence encoding the Influenza virus B (B/PHUKET/3073/2013) hemagglutinin (translated amino acids of EPI529345) (Met1-Thr547), termed as HA, was expressed with a polyhistidine tag at the C-terminus.

### Organism

Influenza B

### Expression Host

Baculovirus-Insect Cells

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE.

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Asp 16

### Molecular Mass

The recombinant hemagglutinin of Influenza virus B (B/PHUKET/3073/2013) consists of 543 amino

acids and predicts a molecular mass of 59.1 kDa.

### Formulation

Lyophilized from sterile 20 mM PB, pH 7.0, 500 mM NaCl, 10 % glycerol.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Influenza B (B/PHUKET/3073/2013)  
Hemagglutinin / HA Protein (His Tag)\n