

SLD5 antibody



Catalog Number: 115355

Product name

SLD5 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

SLD5 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, mouse thymus tissue. Observed molecular weight by Western-blot: 32kd

Preparation

This antibody was obtained by immunization of SLD5 recombinant protein (Accession Number: BC027454). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:115355(SLD5 antibody) at dilution of 1:400