

ITM2C antibody



Catalog Number: 107355

Product name

ITM2C antibody

Specificity

Human; other species not tested.

Antibody description

ITM2C Mouse Monoclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in A431 cells. Observed molecular weight by Western-blot: 35kd

Preparation

This antibody was obtained by immunization of ITM2C recombinant protein (Accession Number: NM_001012516). Purification method: Caprylic acid/ammonium sulfate precipitation.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:

WB: 1:200-1:2000

IF: 1:20-1:200

Validations



Immunofluorescent analysis of HepG2 cells using Catalog No:107355(ITM2C Antibody) at dilution of 1:50 and Rhodamine-Goat anti-Mouse IgG