

# Human TSPAN1 (aa 110-211, Fc Tag) recombinant protein



Catalog Number: 502755

## General Information

### Gene Name Synonym

Tetraspan NET-1; Tetraspanin TM4-C

### Protein Construction

A DNA sequence encoding the human TSPAN1 (O60635)(Tyr110-Asn211) was expressed with the Fc region of rabbit IgG at the N-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

Ser

### Molecular Mass

The recombinant human TSPAN1 comprises 351

amino acids and has a predicted molecular mass of 39.1 kDa. The apparent molecular mass of the protein is approximately 43-53 kDa in SDS-PAGE under reducing conditions due to glycosylation.

### Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human TSPAN1 Protein (aa 110-211, Fc Tag)  
SDS-PAGE