

MPPED2 antibody



Catalog Number: 112755

Product name

MPPED2 antibody

IHC: 1:20-1:200

IF: 1:10-1:100


Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

MPPED2 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells, Hela cells. Positive IHC detected in human brain tissue. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 33kd




Preparation

This antibody was obtained by immunization of MPPED2 recombinant protein (Accession Number: NM_001377954). Purification method: Antigen affinity purified.

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112755(MPPED2 antibody) at dilution of 1:800

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal


Immunofluorescent analysis of HepG2 cells, using MPPED2 antibody Catalog No:112755 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF, IHC



Dilutions

Recommended Dilution:


WB: 1:500-1:5000

Immunohistochemistry of paraffin-embedded

MPPED2 antibody

Catalog Number: 112755

human brain tissue slide using Catalog No:112755(MPPED2 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:112755(MPPED2 Antibody) at dilution of 1:50 (under 40x lens)