

RNF25 antibody

Catalog Number: 114755

Product name

RNF25 antibody

Immunogen

[Human RNF25 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

RNF25 Rabbit Polyclonal antibody. Positive IHC detected in human cervical cancer tissue. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 51 kDa

Preparation

This antibody was obtained by immunization of RNF25 recombinant protein (Accession Number: NM_022453). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114755(RNF25 Antibody) at dilution of 1:3000



Immunohistochemistry of paraffin-embedded human cervical cancer tissue slide using Catalog No:114755(RNF25 Antibody) at dilution of 1:50 (under 40x lens)