

Anti-ZAR1L antibody



Catalog Number: 175755

Product name

Anti-ZAR1L antibody

WB:1:500-2000

IHC-P:1:400-800

Specificity

Human


Validations

Antibody description

Rabbit polyclonal antibody to ZAR1L

Preparation

This antigen of this antibody was klh conjugated synthetic peptide derived from human zar1l 51-150/321




Formulation

Liquid, 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

Fig1: Sample:A549 (Human)cell Lysate at 40 ug; Primary: Anti-ZAR1L at 1/300 dilution; Secondary: IRDye800CW Goat Anti-RabbitIgG at 1/20000 dilution; Predicted band size: 36kD; Observed band size: 36kD

Storage

Store at -20°C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4°C.



Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

IHC-P, WB

Dilutions

Fig2: Paraformaldehyde-fixed, paraffin embedded (Human glioma); Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15min; Block endogenous peroxidase by 3% hydrogen peroxide for 20 minutes; Blocking buffer (normal goat serum) at 37°C for 30min; Antibody incubation with (ZAR1L) Polyclonal Antibody, Unconjugated at 1:400 overnight at 4°C, followed by operating according to SP Kit(Rabbit) (sp-0023) instructions and DAB staining.