

# Anti-Annexin VII Rabbit antibody



Catalog Number: 166755

## Product name

Anti-Annexin VII Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Annexin VII Rabbit polyclonal antibody

## Preparation

Antigen: Recombinant protein of human Annexin VII

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, FC


## Dilutions

WB: 1/2000


IHC: 1/50

FC: 1/20

## Validations



Immunohistochemistry of Annexin VII in paraffin-embedded Human tonsil using Annexin VII Rabbit



Western blot detection of Annexin VII in K562,Rat Brain cell lysates using Annexin VII Rabbit pAb(1:1000 diluted).Predicted band size:53kDa.Observed band size:53kDa.