

# Vinculin antibody



Catalog Number: 116755

---

## Product name

Vinculin antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

Vinculin Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HepG2 cells, NIH/3T3 cells. Observed molecular weight by Western-blot: 117 kDa

## Preparation

This antibody was obtained by immunization of Vinculin recombinant protein (Accession Number: NM\_003373). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116755(Vinculin Antibody) at dilution of 1:2000