

ACVRL1 antibody



Catalog Number: 107755

Product name

ACVRL1 antibody

IP: 1:200-1:1000

IHC: 1:20-1:200

Specificity

Human; other species not tested.

Antibody description

ACVRL1 Rabbit Polyclonal antibody. Positive IHC detected in human lymphoma tissue. Positive FC detected in Jurkat cells. Positive IP detected in Jurkat cells. Positive WB detected in Jurkat cells, human placenta tissue. Observed molecular weight by Western-blot: 63 kDa

Preparation

This antibody was obtained by immunization of ACVRL1 recombinant protein (Accession Number: NM_001406487). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, FC, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

Validations




Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107755(ACVRL1 antibody) at dilution of 1:400




1X10⁶ Jurkat cells were stained with 0.2ug ACVRL1 antibody (Catalog No:107755, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1000.

ACVRL1 antibody


Catalog Number: 107755



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:107755(ACVRL1 antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:107755(ACVRL1 antibody) at dilution of 1:50 (under 10x lens)



IP Result of anti-ACVRL1 (IP:Catalog No:107755, 4ug; Detection:Catalog No:107755 1:300) with Jurkat cells lysate 2400ug.