

Anti-ATF4 Rabbit antibody



Catalog Number: 167065

Product name

Anti-ATF4 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

ATF4 Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human ATF-4

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, FC, IP

Dilutions


WB: 1/1000

IHC: 1/20

FC: 1/20

IP: 1/20

Validations



Western blot detection of ATF4 Rabbit pAb in HeLa cell lysates using ATF4 Rabbit pAb(1:1000 diluted).Predicted band size:39KDa.Observed band size:49KDa.