

Anti-Human IgG Rabbit antibody



Catalog Number: 169065

Product name

Anti-Human IgG Rabbit antibody

Specificity

Human

Antibody description

Human IgG Rabbit polyclonal antibody

Formulation

PBS(pH 7.4) containing with 0.03% Proclin300.

Storage

Store at 4°C short term.Avoid freeze / thaw cycle.

Clonality

Polyclonal


Ig Type

Rabbit IgG

Applications

ELISA

Validations



Indirect ELISA assay for anti-Human IgG rabbit pAb.Antigen coating concentration: 4ug/ml.