

MAN2B2 antibody



Catalog Number: 112465

Product name

MAN2B2 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

MAN2B2 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, mouse lung tissue. Observed molecular weight by Western-blot: 135kd

Preparation

This antibody was obtained by immunization of MAN2B2 recombinant protein (Accession Number: NM_015274). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:112465(MAN2B2 antibody) at dilution of 1:300