

Anti-IKB beta Rabbit antibody



Catalog Number: 167465

Product name

Anti-IKB beta Rabbit antibody

Specificity

Human

Antibody description

IKB beta Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human IKB beta

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, IP

Dilutions


WB: 1/1000

IHC: 1/50

ICC/IF: 1/50

IP: 1/20

Validations



Western blot detection of IKB beta in K562,3T3,HeLa cell lysates using IKB beta Rabbit pAb(1:1000 diluted).Predicted band size:38kDa.Observed band size:48kDa.