

Anti-JNK1 (1H8) Mouse antibody



Catalog Number: 168465

Product name

Anti-JNK1 (1H8) Mouse antibody

Specificity

Human, Mouse

Antibody description

JNK1 (1H8) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human MAPK8 expressed in E. Coli.

Formulation

Purified antibody in PBS with 0.05% sodium azide

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

WB, FC


Dilutions

WB: 1/500 - 1/2000


FC: 1/200 - 1/400

ELISA: 1/10000

Validations



Flow cytometric analysis of HeLa cells using MAPK8 mouse mAb (green) and negative control



Western blot analysis of extracts from Raji, HCT116, Jurkat, MCF7 and A549 cell lysates using JNK1 mouse mAb (1:1000 diluted). Predicted band size: 48 kDa. Observed band size: 46, 54 kDa.