

Integrin beta-6 antibody



Catalog Number: 111865

Product name

Integrin beta-6 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Integrin beta-6 Rabbit Polyclonal antibody. Positive IHC detected in human ovary tumor tissue.

Preparation

This antibody was obtained by immunization of Integrin beta-6 recombinant protein (Accession Number: NM_000888). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:111865(ITGB6 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:111865(ITGB6 antibody) at dilution of 1:50 (under 40x lens)