

# H1N1 HA recombinant protein (C-His)



Catalog Number: 514865

## General Information

### Protein Construction

A DNA sequence encoding the Influenza A virus (A/Fort Monmouth/1/1947) (H1N1) hemagglutinin (ABD77807.1) (Met1-Arg344), termed as HA1, was expressed with a polyhistidine tag at the C-terminus.

### Organism

H1N1

### Expression Host

Human Cells

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE.

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Asp 18

### Molecular Mass

The recombinant HA1 subunit of the Influenza A virus (A/Fort Monmouth/1/1947) (H1N1)

hemagglutinin consists of 338 amino acids and predicts a molecular mass of 38.1 kDa.

### Formulation

Lyophilized from sterile PBS, pH 7.4.

1. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Specific concentrations are included in the hardcopy of COA.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

A hardcopy of COA with reconstitution instruction is sent along with the products. Please refer to it for detailed information.

### SDS-PAGE



H1N1 HA Protein 14111