

Anti-ATXN2L antibody



Catalog Number: 175865

Product name

Anti-ATXN2L antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to ATXN2L

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human atxn2l. aa range:571-620

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of 293T cells using Ataxin-2L Polyclonal Antibody diluted at 1: 2000