

# TBC1D22B antibody



Catalog Number: 115865

## Product name

TBC1D22B antibody

## Specificity

Human; other species not tested.

## Antibody description

TBC1D22B Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

## Preparation

This antibody was obtained by immunization of TBC1D22B recombinant protein (Accession Number: NM\_017772). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:115865(TBC1D22B Antibody) at dilution of 1:500