

Anti-SIRT1 Rabbit antibody



Catalog Number: 166865

Product name

Anti-SIRT1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

SIRT1 Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human SIRT1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of SIRT1 in K562,C6,3T3,HeLa cell lysates using SIRT1 Rabbit pAb(1:1000 diluted).Predicted band size:82kDa.Observed band size:120kDa.