

Anti-PRDM5 antibody



Catalog Number: 101165

Product name

Anti-PRDM5 antibody

Specificity

Human PRDM5

Antibody description

Rabbit polyclonal to PRDM5

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the Human PRDM5, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 10-20 µg/ml

Validations



Items	Lanes	A
Sample (whole cell lysate)		HL-60
Sample Volume (µg/lane)		30
Gel		13% SDS-PAGE reducing gel
Recommended Concentration		10-20 µg/ml
Secondary Antibody		Dylight 800 labeled Antibody to Rabbit IgG (1:1), at 1:5000 dilution.
Explanation		Developed using Odyssey imaging system. Predicted band size: 73 kDa Observed band size: 73 kDa Additional bands at: 56 kDa and 22 kDa (We are unsure as to the identity of these extra bands. They are possible non-specific bindings).

PRDM5 Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot