

Human p63 / TP63 (His & GST Tag) recombinant protein



Catalog Number: 500565

General Information

Gene Name Synonym

Chronic ulcerative stomatitis protein; Keratinocyte transcription factor KET; Transformation-related protein 63; Tumor protein p73-like; p40; p51

Protein Construction

A DNA sequence encoding the human TP63 isoform 1 (Q9H3D4-1) (Met 1-Glu 680) was fused with the N-terminal polyhistidine-tagged GST tag at the N-terminus.

Organism

Human

Expression Host

Baculovirus-Insect Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant human TP63/GST chimera consists of 917 amino acids and has a calculated molecular mass of 105 KDa. It migrates as an approximately 125 KDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile 20mM Tris, 500mM NaCl, pH 7.4, 0.3mM DTT, 20% glycerol

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human P63 / TP63 / Tumor protein p63 Protein (His & GST Tag) SDS-PAGE