

# ZNF251 antibody



Catalog Number: 116965

## Product name

ZNF251 antibody

## Immunogen

[Human ZNF251 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

ZNF251 Rabbit Polyclonal antibody. Positive IF detected in SKOV-3 cells. Positive WB detected in mouse testis tissue, HepG2 cells, mouse skeletal muscle tissue, SKOV-3 cells. Observed molecular weight by Western-blot: 70-76 kDa

## Preparation

This antibody was obtained by immunization of ZNF251 recombinant protein (Accession Number: BC112137). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:116965(ZNF251 Antibody) at dilution of 1:600



Immunofluorescent analysis of SKOV-3 cells using Catalog No:116965(ZNF251 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)