

ALDH1L2 antibody



Catalog Number: 107965

Product name

ALDH1L2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ALDH1L2 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas tissue, human heart tissue. Positive WB detected in mouse kidney tissue, HEK-293 cells. Observed molecular weight by Western-blot: 102 kDa, 89 kDa

Preparation

This antibody was obtained by immunization of ALDH1L2 recombinant protein (Accession Number: BC103935). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:107965(ALDH1L2 antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human pancreas slide using Catalog No:107965(ALDH1L2 Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human pancreas slide using Catalog No:107965(ALDH1L2 Antibody) at dilution of 1:50