

# DNMBP antibody



Catalog Number: 109965

## Product name

DNMBP antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

DNMBP Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, mouse kidney tissue, mouse testis tissue. Positive IP detected in mouse testis tissue,. Observed molecular weight by Western-blot: 180kd

## Preparation

This antibody was obtained by immunization of DNMBP recombinant protein (Accession Number: BC041628). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109965(DNMBP antibody) at dilution of 1:500



IP Result of anti-DNMBP (IP:Catalog No:109965, 4ug; Detection:Catalog No:109965 1:1000) with mouse testis tissue lysate 10000ug.