

HBA1-Specific antibody



Catalog Number: 111265

Product name

HBA1-Specific antibody

Specificity

Human; other species not tested.

Antibody description

HBA1-Specific Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, human spleen tissue. Observed molecular weight by Western-blot: 14-16kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_000558). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:111265(HBA1-Specific antibody) at dilution of 1:300