

ZRANB3 antibody

Catalog Number: 117265

Product name

ZRANB3 antibody

Immunogen

[Human ZRANB3 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

ZRANB3 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Positive IF detected in HEK-293 cells. Observed molecular weight by Western-blot: 150 kDa

Preparation

This antibody was obtained by immunization of ZRANB3 recombinant protein (Accession Number: BC064616). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:117265(ZRANB3 Antibody) at dilution of 1:600



Immunofluorescent analysis of HEK-293 cells using Catalog No:117265(ZRANB3 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)