

Anti-MEK2 Rabbit antibody



Catalog Number: 166665

Product name

Anti-MEK2 Rabbit antibody

Specificity

Human, Mouse, Monkey, Rat

Antibody description

MEK2 Rabbit polyclonal antibody

Preparation

Antigen: Peptide sequence around phosphorylation site of threonine 394 (P-G-T(p)-P-T) derived from Human MEK-2.

Formulation

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB:1/1000

Validations



Western blot detection of MEK2 in C6, COS7, 3T3 and HeLa cell lysates using MEK2 rabbit pAb (1:1000 diluted). Predicted band size: 44KDa. Observed band size: 44KDa.