

Anti-UTF1 antibody



Catalog Number: 176665

Product name

Anti-UTF1 antibody

Ig Type

Mouse IgG2b

Specificity

Human, Mouse, Rat

Applications

WB, ICC

Antibody description

Mouse monoclonal antibody to UTF1

Dilutions

WB: 1:500-1:1,000

Preparation

This antigen of this antibody was peptide

ICC: 1:500-1:1,000

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Validations

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Monoclonal




Fig1: Western blot analysis on D3 cell lysates using anti- UTF1 mouse mAb.